



NRG LIGHT

NrG-AURA Stand Guide

Setup and use guide for the integrated adjustable-angle tabletop stand.

NrG-AURA



Revision 1.0 • August 2026

SETUP

Integrated tabletop stand

The NrG-AURA uses a built-in adjustable-angle stand. No separate floor stand is required for normal tabletop use.

- Place the Aura on a stable, level surface with enough depth to support the complete stand footprint.
- Support the panel with one hand while changing the stand angle. Do not allow the panel to drop or swing freely.
- Set the stand to the desired angle and confirm that the support hardware is fully seated and stable before releasing the panel.
- Route the power cable so it cannot pull the panel toward the edge of the surface or create a trip hazard.
- Keep the rear ventilation area clear.

MOVING THE AURA

Move it as a complete unit

- Turn the panel off and unplug it before moving.
- Allow the panel to cool if it has just completed a session.
- Lift from the panel body/approved handhold areas rather than pulling on the stand or power cord.

Safety: Do not place the Aura on an unstable stool, narrow shelf, soft surface, or any surface where the stand cannot sit flat. If the stand or hinge feels loose, stop use and contact NrG Light.

SUPPORT

Questions about setup?

Contact info@nrglight.com. For panel controls, use the shared NrG Light Panel User Manual.